

NEW DESIGN


Every sampler is vacuum tested for reliable performance.



FEATURES

- Depth markings guide insertion depth and aid in catheter placement
- Atraumatic tip reduces discomfort upon insertion
- Firm and flexible cannula prevents kinking
- New handle design creates a better grip during retrieval of the sample

Endometrial sampler is now non-DEHP & PVC Free

 Not made with natural rubber latex.

TECHNICAL DATA

Depth Mark Range.....	5-11 cm	100% Vacuum Tested
Length	24.5cm	French Size
Sterilization Method	ETO	Box Qty.....
		25 ea.

ORDERING INFORMATION

TMI1175..... Endometrial Sampler

Download IFU



Heartland Medical Specialties, Inc.
8347 Melrose Drive
Lenexa, KS 66214
Telephone: 913 541-1890
FAX: 913 541-1897
E-mail: orders@heartlandmedspec.com