

Designed & manufactured to ensure a minimally invasive yet consistently effective procedure.

FEATURES


Smooth tip reduces risk of trauma upon insertion and removal of catheter.

Stable placement sheath decreases the need for a tenaculum.

Soft, concentric balloon secures above the internal os preventing leakage of sterile saline or water-based contrast media.

Shapeable HS catheter is equipped with an embedded stylet in cases when stentosis is present or when the vaginal tract and cervical canal are not aligned.

Made with 20% Radiopaque material allowing catheter to be seen under xray.

 All Thomas Medical HS Catheters are PVC Free and Non-DEHP.
Not made with natural rubber latex.

Download IFU




TECHNICAL DATA

French Size 5Fr, 7Fr
Length 28cm
Maximum Balloon Volume 5Fr, 1.5ml/7Fr, 3ml
Box Qty. 10 ea.
Sterilization Method R (Gamma)

ORDERING INFORMATION

TMI1154. 5Fr Flexible HS Catheter
TMI1155. 7Fr Flexible HS Catheter
TMI1185. 5Fr Shapeable HS Catheter
TMI1187. 7Fr Shapeable HS Catheter
Bulk shipping options available

SPECIALTY CATHETERS



STEEN OPEN-TIP HS CATHETER

The Steen Open-Tip HS catheter utilizes a balloon closer to the distal tip in order to provide a clearer view of the uterus while performing a hysterosalpingography or sonohysterography procedure.

ORDERING INFORMATION

- TMI1188..... 5Fr Steen Open-Tip HS Catheter
- TMI1189..... 7Fr Steen Open-Tip HS Catheter




MILLER ADVANCE HS CATHETER

The Miller Advance Catheter was designed with features for improved patient comfort with a soft catheter tip and slender catheter shaft that reduces the need for dilation.

ORDERING INFORMATION

- TMI2000..... 5Fr, 21cm Miller Advance HS Catheter
- TMI2001..... 5Fr, 25cm Miller Advance HS Catheter
- TMI2002..... 5Fr, 28cm Miller Advance HS Catheter



TWO-TAP HS CATHETER

The Two-Tap HS catheter utilizes a second stopcock instead of a pinch clamp for easier manipulation when releasing the saline or water-based contrast media from the uterus.

ORDERING INFORMATION

- TMI1181..... 5Fr Two-Tap HS Catheter
- TMI1182..... 7Fr Two-Tap HS Catheter
- TMI1183..... 5Fr Shapeable Two-Tap HS Catheter
- TMI1184..... 7Fr Shapeable Two-Tap HS Catheter

TECHNICAL DATA

CATHETER TYPE	FRENCH SIZE(S)	LENGTH (CM)	MAXIMUM BALLOON VOLUME (ml)	BOX QTY	STERILIZATION METHOD
FLEXIBLE HS CATHETER	5FR/7FR	28	5Fr, 1.5ml/7Fr, 3ml	10	R (GAMMA)
SHAPEABLE HS CATHETER	5FR/7FR	28	5Fr, 1.5ml/7Fr, 3ml	10	R (GAMMA)
STEEN OPEN-TIP HS CATHETER	5FR/7FR	28	5Fr, 1.5ml/7Fr, 3ml	10	R (GAMMA)
MILLER ADVANCE HS CATHETER	5FR	21, 25, 28	5Fr, 1.5ml	10	R (GAMMA)
TWO-TAP HS CATHETER	5FR/7FR	28	5Fr, 1.5ml/7Fr, 3ml	10	R (GAMMA)



Heartland Medical Specialties, Inc.
 8347 Melrose Drive
 Lenexa, KS 66214
 Telephone: 913 541-1890
 FAX: 913 541-1897
 E-mail: orders@heartlandmedspec.com

