


Tested to ensure non-pyrogenicity (LAL) and sperm compatibility



BENEFITS

Shapeability

Shapeable IUI catheters are manufactured with superior memory enabling the doctor to form the tip of the catheter to align with the cervical canal

Clear shapeable IUI catheters allow for full visibility of the washed specimen as it is injected into the uterus

Directional Injection

Side hole near the tip is aligned with the printed depth marks and a flat side of the connector. The alignment allows the catheter to point towards either the left or right fallopian tube ostium for directional injection.

Download IFU



TECHNICAL DATA

French Size 5Fr	Dead Space Volume. see ordering information (ml)
Luer Type Female	Box Qty. 25 ea.

 Not made with natural rubber latex. PVC Free & non-DEHP



ORDERING INFORMATION

Deadspace volumes listed below (.XXml)

TMI1192.....	15cm, 5Fr IUI Catheter (.25ml)
TMI1193.....	15cm, 5Fr Shapeable IUI Catheter (.15ml)
TMI1194.....	18cm, 5Fr IUI Catheter (.28ml)
TMI1195.....	18cm, 5Fr Shapeable IUI Catheter (.20ml)
TMI1295.....	18cm, 5Fr Clear Shapeable IUI Catheter (.20ml)
TMI1196.....	23cm, 5Fr IUI Catheter (.32ml)
TMI1197.....	23cm, 5Fr Shapeable IUI Catheter (.25ml)
TMI1297.....	23cm, 5Fr Clear Shapeable IUI Catheter (.25ml)



Heartland Medical Specialties, Inc.

8347 Melrose Drive

Lenexa, KS 66214

Telephone: 913 541-1890

FAX: 913 541-1897

E-mail: orders@heartlandmedspec.com