



Enhanced features to support visualization of the uterus during pelvic laparoscopic procedures.

IDEAL FOR

Laparoscopic Procedures, Laparoscopic Tubal Occlusion, Diagnostic Laparoscopy, Salpingoplasty

BENEFITS

Reduced catheter shaft diameter (4.5mm) allows for easy insertion

Unique pear-shaped balloon conforms to anatomy and allows for a secure seal of the cervical opening to prevent injection leaks

Flexible extension prevents cracking of luer fitting

Atraumatic beveled catheter tip reduces trauma to the uterine wall.



UNIQUE HANDLE CHARACTERISTICS



TMI6003

Pinch-grip handle allows for quick and effective repositioning during procedures



TMI1151

Snap-on handle allows to easily place the cervical stop against the cervical wall



TMI1171

Sliding handle adjusts after insertion and locks into place to maintain a secure atraumatic hold on the uterus

Handle features larger grip for better leverage during manipulation

ORDERING INFORMATION

TMI1151 Uterine Manipulator/Injector
 TMI6003 K-UMI Manipulator/Injector
 TMI1171 EZ Glide Manipulator/Injector

 Not made with natural rubber latex. PVC Free & non-DEHP

TECHNICAL DATA

UMI O.D./French Size 4.5mm/13.5Fr
 K-UMI O.D./French Size 5.0mm/15Fr
 EZ Glide O.D./French Size 5.0mm/15Fr

Sterilization Method ETO
 Maximum Balloon Volume 10cc
 Length 33cm/13"
 Box Qty. 12 ea.

COMPLEMENTARY PRODUCTS

OS-01-025 os Finders, Box of 25
 US-01-025P Printed Uterine Sound, Box of 25
 Electrodes Various Sizes and Shapes



Heartland Medical Specialties, Inc.

8347 Melrose Drive

Lenexa, KS 66214

Telephone: 913 541-1890

FAX: 913 541-1897

E-mail: orders@heartlandmedspec.com