

ER-SPEC[®]

Single-Use Vaginal Speculum with Built-In Light Source



The medical industry's premier single-use vaginal speculum with built-in light source

OBP Medical's ER-SPEC is a single-use vaginal speculum with built-in LED light source. The ready-to-use device eliminates the time and expense of reprocessing and reduces the risk of cross-contamination in the medical setting by obsoleting rechargeable and plug-in light sources.

ER-SPEC is designed for use on women undergoing a procedure requiring vaginal access and exposure, such as a pelvic examination. Ideal for clinics, emergency rooms, urgent care centers, and physician offices.

OBP
MEDICAL
Simple. Bright. Ideas.



Heartland Medical Specialties, Inc.

Dennis Black
President/Owner

Direct 816-225-9543
Fax 913-541-1897
8347 Melrose Drive
Lenexa, KS 66214

dennis@heartlandmedspec.com
orders@heartlandmedspec.com
www.heartlandmedspec.com

www.obpmedical.com

Think your current light source is clean? Think again.

Think your current light source is clean?

Think again. Recent studies demonstrate a clear potential for transmission of pathogens, and highlight a clear risk to both staff and patients.

- There are 4.5 hospital infections for every 100 patient admissions and nearly 100,000 deaths annually. (CDC)
- There are approximately 1.7 million hospital-acquired infections in U.S. hospitals, costing over \$5 billion annually. (CDC)
- There are bacteria on 86% of instrument handles that are disinfected and designated 'ready for use' in hospital settings. (Journal of Hospital Infection)

TO ORDER:

1-978-291-6853

sales@obpmed.com

www.obpmedical.com

- Item #C020100-1 (Small)
- Item #C020110-1 (Medium)
- Item #C020120-1 (Large)
- Item #C020300 (Small-Sterile)
- Item #C020310 (Medium-Sterile)
- Item #C020320 (Large-Sterile)
- Batteries included in each unit

Product Features:

- Single-use device eliminates the risk of cross-contamination
- Each ER-SPEC comes ready for use with a built-in ultra-bright LED light source
- No assembly required and no additional parts needed
- High powered, battery-operated LED light source can produce peak light for 30+ minutes
- LED light source generates no damaging heat
- Locking hinges offer added strength and stability
- Available in both sterile and non-sterile models



Also available: OfficeSPEC
Single-Use Side-Opening
Vaginal Speculum w/
Built-In LED Light Source.

About Us:

OBP Medical is a leading global developer of single-use, self-contained, illuminating medical devices. Our single-use lights are amongst the brightest on the market, allowing for optimal visualization during procedures. Our innovative products are used in physician offices, surgery centers and more than 3,500 hospitals throughout the U.S., as well as healthcare facilities worldwide.

Medassets Contract # MS02268

HealthTrust Contract # HPG-5741

Premier Contract # PP-AC-120

Amerinet Contract # VH11348

Novation Contract # MS3250

NSN: 6515-01-655-6826

U.S. Patent # 9,307,897

Patent Pending

REV 05282016

OBP Medical Corporation
360 Merrimack Street
Lawrence, MA 01843
www.obpmedical.com

OBP[®]
MEDICAL
Simple. Bright. Ideas.