

# ONETRAC™

Single-Use Cordless Retractor with Built-In Light Source



The medical industry's first cost-effective cordless single-use retractor with built-in light source

Visualization, retraction and hemostasis during a surgical procedure is critical in reducing complications and infections. Avoid the costly and complicated reprocessing requirements of fragile fiber-optic cables. With ONETRAC, access and illumination of the entire pocket is easily and safely achieved.

ONETRAC is an immediate and cost-effective solution for surgical procedures requiring access, illumination and exposure of a surgical tissue pocket or cavity. Illumination in the depths of a surgical tissue pocket or cavity can be met within seconds – no assembly, additional parts/ components, or fiber-optic cables required.

**OBP**  
MEDICAL  
Simple. Bright. Ideas.



*Heartland Medical Specialties, Inc.*

Dennis Black  
President/Owner

Direct 816-225-9543  
Fax 913-541-1897  
8347 Melrose Drive  
Lenexa, KS 66214

dennis@heartlandmedspec.com  
orders@heartlandmedspec.com  
www.heartlandmedspec.com

[www.obpmedical.com](http://www.obpmedical.com)



# Optimal access & visualization without the hassle.

## Tired of expensive retractor illumination systems?

Today's retractor illumination systems cost thousands of dollars with hundreds more in ongoing disposable and/or reprocessing costs. ONETRAC is the first line of cost-effective cordless single-use surgical retractors complete with built-in LED light source.

Each ONETRAC comes ready for use with everything needed to successfully retract soft tissue and illuminate a surgical pocket or cavity all at a fraction of the cost of other retractor illumination systems. No assembly, additional parts/components, or fiber-optic cables are required.

## TO ORDER:

1-978-291-6853

sales@obpmed.com

www.obpmedical.com

- Item # C100110 (135mm x 30mm)
- Item # C100140 (90mm x 22mm)  
(Available September 2016)
- Item # C100150 (40mm x 20mm)  
(Available October 2016)
- Batteries included in each unit
- Each unit comes packed sterile

## Product Features:

- Disposable, single-use device eliminates risk of cross-contamination
- On/off switch to preserve and maximize duration of peak illumination strength
- Built-in LED light source produces no damaging heat and will provide peak power for 100 + minutes
- Built-in smoke evacuation channel
- Revolutionary reinforced nylon material provides exceptional strength and rigidity
- Angled tip with blunt teeth for improved tissue grip
- Lightweight for added comfort and ease of use



### About Us:

OBP Medical is a leading global developer of single-use, self-contained, illuminating medical devices. Our single-use lights are amongst the brightest on the market, allowing for optimal visualization during procedures. Our innovative products are used in physician offices, surgery centers and more than 3,500 hospitals throughout the U.S., as well as healthcare facilities worldwide.

OBP Medical Corporation  
360 Merrimack Street  
Lawrence, MA 01843  
www.obpmedical.com

**OBP**<sup>®</sup>  
MEDICAL  
Simple. Bright. Ideas.