


• COMEG scopes and imaging products are built with superior quality, commitment, and dedication.



- Rod lense, objective and ocular are multilayer **anti-reflective coated** for increased optical performance



- **Fiber illumination** for Field of View optimization



- **Three tube** design provides increased strength and durability

- **Laser welded** manufacturing ensures precision and quality

- Proximal and distal window made of **sapphire** to ensure for scratch resistance

- **Gold soldering** of the proximal and distal end provides increased durability and lifespan throughout continued sterilization

- **Stainless steel** tubes and scope body allow for optimal quality and integrity

## Discover our additional Equipment Platforms



**SYMBIOZ**  
FULL HD multidisciplinary endoscopy tower



**ULYSS**  
Multi-view Sinuscope



**UBIPACK**  
Mobile office endoscopy solution for GYN, URO and ENT



**PEZOTOME M+**  
Piezo-Ultrasonic Surgery Unit  
Ultrasonic piezo bone surgery platform

## • Our vision

<b>ALWAYS...</b>	<b>EVERYWHERE...</b>	<b>AT YOUR SIDE.</b>
<b>Dedicated</b> to medical advancements for more than 40 years.	<b>Worldwide</b> centers allow for immediate product availability and rapid service.	<b>Efficient</b> communication with our local sales force, and nationally placed inventory, to ensure reduced delivery times.

## • 40 years of expertise

- 1975** COMEG founded in Tuttlingen, Germany
- 1977** SOPRO founded in La Ciotat, France
- 1977** Introduction of medical endoscopes
- 1986** Micro-cameras designed specifically for endoscopy
- 2001** Engineered and introduced a new generation of high quality medical optics
- 2006** COMEG and SOPRO merged
- 2009** Introduction of the SOPRO COMEG Full HD series of cameras combined with SOPRO Medical Software
- 2011** Introduction of ULYSS multi-view sinuscope and UBIPACK; new best in class solutions for ENT and GYN
- 2015-2016** International launch of SYMBIOZ; a multidisciplinary, HD, surgical endoscopy tower  
International launch of M+ ultrasonic piezo bone surgery platform

## • Production sites

### La Ciotat - France



Medical Imaging

### Tuttlingen - Germany



Endoscopes and Accessories

### Mérignac - France



M+

**HEAD OFFICE**  
SOPRO SA  
ZAC Athélia IV  
Avenue des Gènévriers  
13705 La Ciotat cedex  
France  
Tel. +33 [0] 442 98 01 01  
info@comegmedical.com

**GERMAN OFFICE**  
SOPRO-COMEG GmbH  
Dornierstr. 55  
D-78532 Tuttlingen  
Germany  
Tel. +49 [0] 7461 9344 0

**U.S.A. OFFICES**  
124 Gaither Drive Suite 140  
Mount Laurel, NJ 08054  
U.S.A.  
Tel. +1 856 222 9988


6600 France Avenue South  
Edina, Mn. 55435  
U.S.A.  
Tel. +1 844 99COMEG

www.comegmedical.com  
orders.usa@comegmedical.com


**COMEG**  
medical technologies

Reference document: Non contractual document - ref. 011726A - 02/2017 - Copyright © 2017 SOPRO. All rights reserved. COMEG Medical Technologies is a Division of ACTEON Group. No information or part of this document may be reproduced or transmitted in any form without the prior permission of SOPRO.

# Over 40 years of expertise in the design and production of endoscopes



**COMEG**  
medical technologies




# GYNECOLOGY

Clarity and efficiency are ideals we strive for with our range of Gynecology products.



Product Code	Description
<b>HYSTEROSCOPES</b>	
163 282 990	2.7mm, 0° Wide Angle Hysteroscope, 301mm length
163 282 992	2.7mm, 30° Wide Angle Hysteroscope, 301mm length
165 293 000	2.9mm, 0° Hysteroscope, 301mm length
165 293 001	2.9mm, 12° Hysteroscope, 301mm length
165 293 002	2.9mm, 30° Hysteroscope, 301mm length
163 412 990	4mm, 0° Wide Angle Hysteroscope, 301mm length
163 412 981	4mm, 12° Hysteroscope, 301mm length
163 412 992	4mm, 30° Wide Angle Hysteroscope, 301mm length



# LAPAROSCOPY

Our Full HD line of Laparoscopes provides superior image quality while remaining an extremely cost effective solution.



Product Code	Description
<b>LAPAROSCOPES</b>	
162 053 150	5.5mm, 0° Full HD Laparoscope, 316mm length
162 053 152	5.5mm, 30° Full HD Laparoscope, 316mm length
162 053 153	5.5mm, 45° Full HD Laparoscope, 316mm length
162 103 250	10mm, 0° Full HD Laparoscope, 323mm length
162 103 252	10mm, 30° Full HD Laparoscope, 325mm length
162 103 253	10mm, 45° Full HD Laparoscope, 325mm length



# UROLOGY

At COMEG we provide broad product offerings for the most demanding Urology procedures.



Product Code	Description
<b>CYSTOSCOPES</b>	
165 293 000	2.9mm, 0° Cystoscope, 301mm length
165 293 001	2.9mm, 12° Cystoscope, 301mm length
165 293 002	2.9mm, 30° Cystoscope, 301mm length
165 403 000	4mm, 0° Cystoscope, 301mm length
165 403 001	4mm, 12° Cystoscope, 301mm length
165 403 002	4mm, 30° Cystoscope, 301mm length
165 403 005	4mm, 70° Cystoscope, 303mm length
<b>URETERO-RENOSCOPES</b>	
186 552 001	Uretero-renaloscope 7.5-6 Fr., 12° Bichler style, 315mm length instrument Channel 30.000 pixels
186 528 001	Uretero-renaloscope 7.5-6 Fr., 12° Bichler style, 430mm length instrument Channel 30.000 pixels
186 395 001	Uretero-renaloscope 8.9-8 Fr., 12° Bichler style, 315mm length instrument Channel 30.000 pixels
186 403 001	Uretero-renaloscope 8.9-8 Fr., 12° Bichler style, 430mm length instrument Channel 30.000 pixels



# ENT

Whether in an office or operating room, COMEG offers a wide variety of ENT products to meet various clinical needs.



Product Code	Description
<b>SINUSCOPES</b>	
161 271 780	2.7mm, 0° Wide Angle Sinuscope, 175mm length
161 271 782	2.7mm, 30° Wide Angle Sinuscope, 175mm length
161 271 785	2.7mm, 70° Wide Angle Sinuscope, 175mm length
161 401 810	4mm, 0° Wide Angle Sinuscope, 175mm length
161 401 812	4mm, 30° Wide Angle Sinuscope, 175mm length
161 401 813	4mm, 45° Wide Angle Sinuscope, 175mm length
161 401 815	4mm, 70° Wide Angle Sinuscope, 175mm length
<b>OTOSCOPES</b>	
161 300 650	3mm, 0° Otoscope, 69mm length
161 400 650	4mm, 0° Otoscope, 69mm length
<b>LARYNGOSCOPES</b>	
162 081 685	8.5mm, 70° Laryngoscope, 185mm length
162 081 687	8.5mm, 90° Laryngoscope, 185mm length
<b>ULYSS</b>	
161 401 760	4.3mm, 10° to 100° view FULL HD Multi-view Sinuscope, 0° to 340° shaft rotation, titanium body, 175mm length

Discover additional COMEG product offerings for



Arthroscopy



Spine




CMF



Capital Equipment

Our vision is our value. *Always...Everywhere...At your side.*



- 4.3mm
- 10° to 100°
- 340° rotation
- 175mm working length
- Autoclavable
- Titanium body

Request our catalog  
Contact your local representative

Visit our website

[www.comegmedical.com](http://www.comegmedical.com)